

# Bedford Wrist/Thumb Brace

A unique wrist immobilisation brace which can be used to support un-displaced fractures of the scaphoid, trapezium and thumb metacarpal. It features a mouldable bar which holds the thumb firmly in abduction and has been improved to reduce pressure points. It also now features an open thumb seam for ease of application and comfort.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Instability
- + Post fracture management





Immobilises wrist and holds thumb firmly in abduction with a mouldable metal bar.



The Bedford Wrist/Thumb Brace also now features an open thumb seam for ease of application and comfort.



Can be used to support undisplaced fractures of the scaphoid, trapezium and thumb metacarpal.

Code (Left)	Code (Right)	Size
BWTB/12	BWTB/22	XS
BWTB/13	BWTB/23	S
BWTB/15	BWTB/25	M
BWTB/17	BWTB/27	L
BWTB/18	BWTB/28	XL