

# Single Armchair Finger Splint

A dynamic splint with coil design for PIP and DIP joint flexion contractures and Boutonnieres deformity. Maintains assisted extension at the PIP joint while allowing active flexion and can be used for soft tissue injury affecting the PIP joint with a possible flexion deformity.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + PIP joint flexion contractures
- + DIP joint flexion contractures
- + Boutonniere deformity
- + Soft tissue injury affecting the PIP joint






The dynamic spring supports the PIP joint, leaving the DIP free to flex.



Maintains assisted extension at the PIP joint and allows active flexion of the finger.



Can be used for soft tissue injury affecting the PIP joint. Also ideal for Boutonniere and other flexion deformities.

Code	Size
CAP/2XS	XS
CAP/2S	S
CAP/2M	M
CAP/2L	L
CAP/2XL	XL