

Deluxe Maternity Belt

Anterior foam laminate panels - low profile shape with hook / loop closure for improved fit, especially when sitting.

Adjustable stretch fulcrums give additional support. The Deluxe Maternity Belt helps relieve lower back pain and associated symptoms.



Code	Circumference (at widest part of abdomen)
DMB/1	90 - 96cm
DMB/2	96 - 102cm
DMB/3	102 - 108cm
DMB/4	108 - 116cm
DMB/5	116 - 126cm
DMB/6	126 - 136cm
DMB/7	136 - 146cm
DMB/8	146 - 156cm
DMB/9	156 - 165cm
DMB/10	165 - 178cm