

# Duofit Maternity Belt

20cm deep circumferential elastic wrap with anterior cut-away to improve comfort when sitting. Angled hook/loop closure to allow accurate positioning under tummy.

Additional narrow pelvic wrap provides extra support around the pelvis for improved stability at SI joints.

Also available in black.



Code (White)	Code (Black)	Circumference (at widest part of the abdomen)
DUB/1	DUB/1/BK	90 - 96 cm
DUB/2	DUB/2/BK	96 - 102cm
DUB/3	DUB/3/BK	102 - 108cm
DUB/4	DUB/4/BK	108 - 116cm
DUB/5	DUB/5/BK	116 - 126cm
DUB/6	DUB/6/BK	126 - 136cm
DUB/7	DUB/7/BK	136 - 146cm
DUB/8	DUB/8/BK	146 - 156cm
DUB/9	DUB/9/BK	156 - 165cm
DUB/10	DUB/10/BK	165 - 179cm