

Paediatric High Sided Cast Sandal

Designed with a hard-wearing non-stretch nylon-faced upper for long-term wear, the Paediatric High Sided Cast Sandal accommodates bulky casts and dressings to prevent any delay in healing. The sandal features a loop lock D-ring fastening for secure and easy fastening, as well as a non-slip rocker sole for stability.



Code	Sole Length x Width	Size
HSCSL/CS	14cm x 5.5cm	S
HSCSL/CM	17cm x 6.5cm	M
HSCSL/CL	20cm x 7.5cm	L