

Jura Black Wrist Brace

The Jura Black Wrist Brace immobilises and protects the wrist following tendonitis and other soft tissue injuries such as sprains/ strains, and post-trauma healing. Lightweight, breathable and contoured for a comfortable fit. Slide-on design with D-ring closures allow easy application.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Stable fractures
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Arthritis
- + Peripheral nerve pathologie.
- + Slide-on design with D-ring closures allow easy application especially for patients with reduced dexterity.





Made from dual compression elastic for a great fit, the Neoprene fabric is breathable and comfortable against the skin.



Durable and tough webbing bar casing and a two-piece brace with a longitudinal seam conforms perfectly to the wrist.



Hook and loop fastenings make it easy to apply and remove.

Code (Left)	Code (Right)	Size
JRWB/12	JRWB/22	XS
JRWB/13	JRWB/23	S
JRWB/15	JRWB/25	M
JRWB/17	JRWB/27	L
JRWB/18	JRWB/28	XL