

# Jura Black Wrist Brace - Long

Similar design and function as the standard Jura Wrist Brace. Increasing length to 23cm offers improved wrist stability and immobilisation. Dual compression elastic for improved support and fit. A two-piece wrist brace with longitudinal seam conforms perfectly to the wrist, hand and forearm. Aluminium palmar bar and thumb stay provides excellent support to wrist and first CMC & MCP joints. The thumb stay can be easily adjusted for the required degree of abduction and the thumb closure is easily adjusted for a comfortable fit.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Stable fractures
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Arthritis
- + Peripheral nerve pathologies





Made from dual-compression elastic with longitudinal seams that conform to the wrist, hand and forearm.



Slide-on design with D-ring closures allow for easy application, especially for patients with reduced dexterity.

Code (Left)	Code (Right)	Size
JRWBL/12	JRWBL/22	XS
JRWBL/13	JRWBL/23	S
JRWBL/15	JRWBL/25	M
JRWBL/17	JRWBL/27	L
JRWBL/18	JRWBL/28	XL