

Jura

Wrist Thumb

Dual compression elastic for improved support and fit. A two-piece wrist thumb brace with longitudinal seam conforms perfectly to the wrist, hand and forearm. Aluminium palmar bar and thumb stay provides excellent support to wrist and first CMC & MCP joints. The thumb stay can be easily adjusted for the required degree of abduction and the thumb closure is easily adjusted for a comfortable fit.



Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Carpal tunnel syndrome
- + CMC instability
- + De Quervain's syndrome
- + Gamekeeper's thumb
- + General weakness
- + Guyon's canal syndrome
- + Inflamed/swollen tendons
- + Luxation
- + Overuse injuries
- + Post carpal/metacarpal surgery
- + Post surgical repair
- + Post traumatic conditions
- + Repetitive strain injury
- + Rhizarthrosis
- + Scapholunate instability
- + Scapholunate tear
- + Skier's thumb
- + Soft tissue injury
- + Sprained wrist + more



Code (Left)	Code (Right)	Size
JWT/12	JWT/22	XS
JWT/13	JWT/23	S
JWT/15	JWT/25	M
JWT/17	JWT/27	L
JWT/18	JWT/28	XL