

Junior Wrist/Thumb Brace

A wrist brace made of Neoprene which can be cut without fraying to fit around swollen joints, giving compression and excellent conformity.

Incorporating a malleable metal thumb stay to control an unstable CMC joint; the thumb section can be trimmed to allow IP flexion.


Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Juvenile chronic arthritis
- + Soft tissue injury





Can be trimmed without fraying to fit around swollen joints, and provides compression and conformity to ease pain and discomfort.



Malleable thumb stay helps to control an unstable CMC joint, thumb section can also be trimmed to allow IP flexion.



It provides compression and conforms well around the joints, easing pain and discomfort while the area heals.

Code (Left)	Code (Right)	Size
JWTO/13	JWTO/23	S
JWTO/15	JWTO/25	M
JWTO/17	JWTO/27	L