

# LMB Air-soft Resting Hand Splint

Maintains mobility while providing support for the hand. The malleable aluminum frame with a special self-wicking cover helps control perspiration dampness.

Supports the fingers, thumb and wrist in a functional resting position while remaining easily adjustable. Allows up to 60° MP flexion and wrist flexion/extension. Indicated for use in functional positioning for arthritis, post-surgery or stroke.

To size, measure width at the MP joints.



Code	Style	Size
L/325-BL	Left • 5.7cm - 7cm	S
L/325-BR	Right • 5.7cm - 7cm	S
L/325-CL	Left • 7.3cm - 8.3cm	M
L/325-CR	Right • 7.3cm - 8.3cm	M
L/325-DL	Left • 8.6cm - 9.5cm	L
L/325-DR	Right • 8.6cm - 9.5cm	L