

Lumbarform Belt with Stays

Triple elastic fulcrum straps at the back which attach to one simple strap at the front, giving extra compression on the spinal curve.

Two flexible posterior stays and a double-layer ventilated elastic back pocket stabilise the spine. The ventilated elastic which has Lycra binding to prevent chafing.

Lumbarform Belt with Stays has a back depth of 23cm.



Code	Size	Waist Measurement
LBS/93	S	69 - 76 cm
LBS/95	M	79 - 86 cm
LBS/97	L	89 - 97 cm
LBS/98	XL	99 - 107 cm