

Liberty™ Flare Wrist Brace

The Liberty™ Flare Wrist Brace is designed specifically for people whose wrist circumference and mid-forearm circumference vary by 10cm or more. Made of woven cotton and coated rubber. Ideal for the individual whose wrist and hand are small in proportion to the forearm.

One adjustable size fits most large forearms. Wrist circumference 15cm - 22cm, forearm circumference 25cm - 36cm.



Code	Size
NC/14562	Left • Regular
NC/14563	Right • Regular
NC/14560	Left • Large
NC/14561	Right • Large