

# Liberty™ Elastic Wrist Brace

The Liberty™ Elastic Wrist Brace allows the fingers and thumb to move freely while providing excellent wrist support. The unique cutout in the aluminium stay allows for unimpeded range of motion of the thumb. Anatomically designed stay can be modified easily for a custom fit. Made of woven cotton and coated rubber, the open-weave stretch fabric allows for air circulation and helps prevent skin maceration.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Inflammation
- + Strains
- + Carpal tunnel syndrome

To size, measure wrist circumference.



Code (Beige)	Wrist Circumference	Size
NC/15871	Left • 11cm - 14cm	XS/Youth
NC/15872	Right • 11cm - 14cm	XS/Youth
NC/15873	Left • 14cm - 17cm	S
NC/15874	Right • 14cm - 17cm	S
NC/15875	Left • 17cm - 18cm	M
NC/15876	Right • 17cm - 18cm	M
NC/15877	Left • 18cm - 20cm	L
NC/15878	Right • 18cm - 20cm	L
NC/15879	Left • 20cm - 22cm	XL
NC/15880	Right • 20cm - 22cm	XL

Code (Blue & Turquoise)	Wrist Circumference	Size
NC/17671	Left • 11cm - 14cm	XS/Youth
NC/17672	Right • 11cm - 14cm	XS/Youth
NC/17673	Left • 14cm - 17cm	S
NC/17674	Right • 14cm - 17cm	S
NC/17675	Left • 17cm - 18cm	M
NC/17676	Right • 17cm - 18cm	M
NC/17677	Left • 18cm - 20cm	L
NC/17678	Right • 18cm - 20cm	L
NC/17679	Left • 20cm - 22cm	XL
NC/17680	Right • 20cm - 22cm	XL

Code (Blue & Lavender)	Wrist Circumference	Size
TWB/CL/13	Left • 11cm - 14cm	XS/Youth
TWB/CL/23	Right • 11cm - 14cm	XS/Youth
TWB/CL/15	Left • 14cm - 17cm	S
TWB/CL/25	Right • 14cm - 17cm	S
TWB/CL/17	Left • 17cm - 18cm	M
TWB/CL/27	Right • 17cm - 18cm	M
TWB/CL/18	Left • 18cm - 20cm	L
TWB/CL/28	Right • 18cm - 20cm	L
TWB/CL/19	Left • 20cm - 22cm	XL
TWB/CL/29	Right • 20cm - 22cm	XL