

Progress[™] Palm Protector

Trim and shape for truly custom progressive positioning. Progress[™] Palm Protector Orthosis is made 9.5mm thick, high-density foam covered with soft fabric and a hook and loop closure.

Removable finger separators offer finger abduction to assist with tone control and hygiene. Fold or roll the palmar flap to achieve the desired thickness.

When rolled, the flap can accommodate a small cone or cylindrical foam to ease the hand out of a flexed posture.

One size fits most. Can be trimmed for a custom fit.



Code	Style
NC/25620	Left
NC/25621	Right
NC/25629	Replacement Finger Separator