

# Comfort Cool™ D-ring Wrist Support (Long)

Cool, lightweight support provides cumulative trauma relief. Made of 1.6mm perforated neoprene and lined with Terry cloth, the Comfort Cool™ D-ring Wrist Support (Long) provides excellent support without excessive bulk or warmth.

The aluminium volar support is pre-bent for 15° wrist extension and can be removed for fitting adjustments. Circumferentially designed with a double-layer cotton stockinette under D-ring closures providing comfort and easy application. Anti-Microbial and Latex free.

Attach Comfort-Gel Pads, for added scar care or padding. To size, measure wrist circumference. 23cm in length.





The aluminium volar support is pre-bent for 15° wrist extension and can be removed for fitting adjustments.



Circumferentially designed with a double-layer cotton stockinette under D-ring closures providing comfort and easy application.

| Code     | Wrist Circumference | Size |
|----------|---------------------|------|
| NC/52970 | Left • < 15cm       | XS   |
| NC/52971 | Right • < 15cm      | XS   |
| NC/52972 | Left • 15cm - 17cm  | S    |
| NC/52973 | Right • 15cm - 17cm | S    |
| NC/52974 | Left • 17cm - 19cm  | M    |
| NC/52975 | Right • 17cm - 19cm | M    |
| NC/52976 | Left • 19cm - 22cm  | L    |
| NC/52977 | Right • 19cm - 22cm | L    |
| NC/52978 | Left • 22cm >       | XL   |
| NC/52979 | Right • 22cm >      | XL   |