

# Comfort Cool<sup>™</sup> Thumb CMC Abduction Splint

Comfortable splint provides direct support for the thumb CMC joint while allowing full finger function. The same CMC thumb wrap-around orthosis as the Comfort Cool<sup>®</sup> CMC Restriction Splint, but with the addition of a firm, moldable insert inside the thumb strap.

Thumb strap can be microwaved to soften the thermoplastic insert, and then molded to fit precisely along the thenar eminence crease. The thumb strap helps to abduct and stabilize the first metacarpal and open the thumb web space.



Code	DPC Circumference	Size
NC/79490	Left • 11cm - 14cm	Toddler
NC/79491	Right • 11cm - 14cm	Toddler
NC/79492	Left • 14cm - 17cm	Youth
NC/79493	Right • 14cm - 17cm	Youth
NC/79494	Left • 15cm - 18cm	Small
NC/79495	Right • 15cm - 18cm	Small
NC/79496	Left • 18cm - 20cm	Medium
NC/79497	Right • 18cm - 20cm	Medium
NC/79498	Left • 20cm - 23cm	Large
NC/79499	Right • 20cm - 23cm	Large