

Comfort Cool™ Thumb Abductor Strap

Made from 3.2mm thick neoprene strap that provides a comfortable stretch for the finger or thumb MP's, designed to treat thumb spasticity, tendinitis, or trigger finger.

Can be used to decrease thumb adductor tone or to position the thumb for optimal hand function. Achieve radial or palmar thumb abduction by relocating the thumb loop on the separate wrist strap.

To size, measure circumference at the wrist crease. Small size is appropriate for some children.



Code	MP Circumference	Size
NC/79520	11cm - 14cm	S
NC/79521	Right • 16cm - 19.5cm	M
NC/79522	Left • 19.5cm - 22.5cm	L