

Comfort Cool™ Thumb Spica Splint

A firm yet flexible thumb and wrist support, the Liberty™ Wrist & Thumb Spica secures thumb IP and MP joints. Made of lightweight trilaminate foam with wide stretch loop straps and soft lining, it offers protection, support and positioning of the thumb and wrist for painful and inflammatory conditions.

It's perforated foam contours to the forearm and provides ventilation for a cool, comfortable wear. Easy modification with scissors and a heat gun. Material will not fray. Can be used for post-surgical and post-fracture positioning. Extended thumb design holds the thumb IP and MP joints securely. Approximately 30cm long. Latex free.



Code	Style	Size
NC/79528	Short	S/M
NC/79529	Short	M/L
NC/79530	Mid-Length	S/M
NC/79531	Mid-Length	M/L
NC/79532	Long	S/M
NC/79533	Long	M/L