

Comfort Cool™ D-ring Thumb & Wrist Splint

Thumb tip can be trimmed to allow full IP flexion. Made of thin, 1.6mm perforated neoprene with terry cloth liner that helps keep skin dry. The neoprene can be trimmed and will not fray. The Comfort Cool™ D-ring Thumb & Wrist Splint protects and supports the entire thumb and wrist for conditions such as arthritis, tendinitis, sprains, strains, skier's thumb, gamekeeper's thumb, De Quervain's, wrist instabilities and post-surgical repair.

Metal volar wrist stay is pre-bent to position the wrist in functional extension (approximately 15°). Volar wrist and thumb stays can be bent to conform as needed and can be removed as the condition improves. Two lengths available: Long measures 5.1cm longer than Regular. Anti-microbial and latex free.





Made of thin, 1.6mm perforated neoprene with terry cloth liner that helps keep skin dry.



Metal volar wrist stay is pre-bent to position the wrist in functional extension (approximately 15°). Volar wrist and thumb stays can be bent to conform as needed and can be removed as the condition improves.

Code (Regular)	Code (Long)	Wrist Circumference	Size
NC/91250	NC/91260	Left • 15cm - 18cm	XS
NC/91251	NC/91261	Right • 15cm - 18cm	XS
NC/91252	NC/91262	Left • 18cm - 20cm	S
NC/91253	NC/91263	Right • 18cm - 20cm	S
NC/91254	NC/91264	Left • 20cm - 21cm	M
NC/91255	NC/91265	Right • 20cm - 21cm	M
NC/91256	NC/91266	Left • 21cm - 23cm	L
NC/91257	NC/91267	Right • 21cm - 23cm	L
NC/91258	NC/91268	Left • 23cm >	XL
NC/91259	NC/91269	Right • 23cm >	XL