

Proform Belt

An elastic front with hook/loop closure and a cotton back section, with detachable comfort pad and two malleable steels which can be shaped to fit the spinal curve.

A webbing strap which can be positioned in the lumbar or sacral area for extra compression on the spinal curve. Back depth is 28cm to support the thoracic area reducing to 20cm at the front for comfort.



Code	Waist Measurement	Size
PFM/83	58 - 66 cm	S
PFM/85	68 - 76 cm	M
PFM/87	79 - 86 cm	L
PFM/88	89 - 97 cm	XL
PFM/80	99 - 107 cm	XXL