

Proflex Belt

Hook/loop fastening elastic belt, with an elastic back pocket to accommodate a panel. Available with a choice of back panels.

A cotton fulcrum strap which can be positioned in the lumbar or sacral area from extra compression on the spine. The front depth is 20cm.

Flexible panel option - The back pocket encloses a panel of polythene lined with closed cell foam, which becomes flexible with body heat.



Code (No Panel)	Code (Flexible Panel)	Size	Waist Measurement
PFN/83	PF/83	S	64 - 76cm
PFN/85	PF/85	M	79 - 89cm
PFN/87	PF/87	L	89 - 99cm
PFN/88	PF/88	XL	99 - 109cm
PFN/89	PF/89	XXL	109 - 119cm