

Poroflex Pull & Thumb

The Poroflex Pull & Thumb Brace is a simple slip-on design with two strap D-ring closures. Excellent for patients with limited dexterity.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Stable fractures
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Arthritis
- + Peripheral nerve pathologies
- + Slide-on design with D-ring closures allow easy application especially for patients with reduced dexterity





Poroflex is a new three dimensional fabric offering improved patient comfort due to breathability and temperature control.



Features a thumb extension for MCP & CMC control as well as an adjustable aluminium thumb stay.



Has an adjustable palmar bar with flexible dorsal stays for excellent wrist compression and stability.

Code (Left)	Code (Right)	Size
PFTB/12	PFTB/22	XS
PFTB/13	PFTB/23	S
PFTB/15	PFTB/25	M
PFTB/17	PFTB/27	L
PFTB/18	PFTB/28	XL