

Paediatric Wrist Splint

Providing adequate support is crucial to the healing process, whether from surgery, arthritic condition or an injury. The Paediatric Lycrafleece Thumb Spica includes a flexible splint to provide thumb stabilisation and is made from comfortable, breathable Lycrafleece material.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis of the thumb
- + Basal joint arthritis
- + CMC arthritis
- + CMC instability
- + Controls thumb abduction
- + Gamekeeper's thumb
- + Inflamed/swollen tendons
- + Juvenile Chronic arthritis
- + Luxation
- + Medial subluxation
- + Osteoarthritis
- + Post carpal/metacarpal surgery
- + Post surgical repair
- + Provides radial/ulnar support
- + Repetitive Strain Injury (RSI)
- + Rheumatoid arthritis
- + Rhizarthrosis
- + Scapholunate instability
- + Skier's thumb
- + Stener Lesion
- + Tendinopathy
- + Tendomyopathy
- + Worn cartilage
- + Wrist Arthralgia
- + Wrist contusion
- + Wrist instability/strain
- + Wrist tendonitis
- + Wrist trauma





Lycrafleece is a Lycra Spandex 3-ply laminate, designed for luxurious comfort and impressive durability. The Lycra Spandex gives the material the ideal amount of stretch to support the thumb without pulling on it while the polyester lining material wicks moisture away from the skin to stop sweat or heat from building up under the thumb brace.



Includes a flexible splint to provide thumb stabilisation, while the outer material is soft, long-lasting, elasticated and extremely breathable.



Code (Left)	Code (Right)	Size
PLWB/12	PLWB/22	XS
PLWB/13	PLWB/23	S
PLWB/15	PLWB/25	M
PLWB/17	PLWB/27	L
PLWB/18	PLWB/28	XL