

Pro-Rheuma Thumb Spica

Made from a breathable high compression material with microfleece lining and metal stay for added rigidity. The Pro-rheuma Thumb Spica features a removable semi-flexible thumb stay that helps stabilise CMC and MCP joints.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Instability





Breathable, high compression material with microfleece lining.



Removable semi-flexible thumb stay helps stabilise CMC & MCP joints.

Code (Left)	Code (Right)	Size
PRHS/12	PRHS/22	XS
PRHS/13	PRHS/23	S
PRHS/15	PRHS/25	M
PRHS/17	PRHS/27	L
PRHS/18	PRHS/28	XL