

Procool

Wrist Thumb Brace

An adjustable wrist thumb wrap offering good support around the wrist, MCP and CMC joints. Adjustable thumb closure with excellent shaping around the palm.

Procool is a range of lightweight and breathable splints producing excellent compression and support for a wide range of conditions and injuries. Made with tiny perforations in the Neoprene for comfortable wear and a terry cloth lining that feels good against the skin. Procool offers you just the fit you want - excellent stability without excessive bulk.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Basal thumb arthritis
- + Thumb CMC & MCP arthritis
- + Post UCL injury
- + Degenerative joints
- + Rheumatoid arthritis
- + Osteo-arthritis





Made with tiny perforations in the Neoprene for comfortable wear and a terry cloth lining that feels good against the skin. It also provides circumferential compression, helping to improve circulation and reduce swelling.



Designed to stabilise the appropriate joints while allowing you to maintain the manual dexterity necessary for so much of day-to-day life.

Code (Left)	Code (Right)	Size
PWTS/12	PWTS/22	XS
PWTS/13	PWTS/23	S
PWTS/15	PWTS/25	M
PWTS/17	PWTS/27	L
PWTS/18	PWTS/28	XL