

SK Knee Brace

The front opening design and easy-open tabs on the hook/loop straps make this brace suitable for the arthritic or elderly patient, who may suffer from limited dexterity.

Hyperstop hinges provide medial and lateral support for the unstable knee, while allowing flexion. Also available with autolock hinge.

- + Mild instabilities of the knee
- + Meniscus injuries
- + Trauma or degenerative deformities of the knee joint
- + Protection in the later stages of rehab



Code	Knee Circumference	Size
SK/02	25cm - 28cm	XS
SK/03	28cm - 34cm	S
SK/05	34cm - 40cm	M
SK/07	40cm - 46cm	L
SK/08	46cm - 56cm	XL
SK/09	56cm - 61cm	XXL