

Stabilised Knee Sleeve

A pull-on Neoprene sleeve incorporating a buttress pad, spiral stays for medial and lateral support, and hook/loop straps above and below the knee for improved fit. The Stabilised Knee Sleeve can be selected with a buttress pad to oppose medial/lateral patella migration; caused by ligament deficiency.

The buttress pad positions are defined as viewed by the patient, for example a left buttress pad would give lateral support to the left leg and medial support to the right leg. A full buttress is also available.

+ Mild instabilities of the knee



Code (Left)	Code (Right)	Circumference 15cm Below Knee	Knee Circumference	Circumference 15cm Above Knee	Size
ST/12	ST/22	30cm	30cm - 33cm	38cm	XS
ST/13	ST/23	33cm	33cm - 36cm	46cm	S
ST/15	ST/25	38cm	36cm - 38cm	50cm	M
ST/17	ST/27	41cm	38cm - 43cm	53cm	L
ST/18	ST/28	43cm	43cm - 48cm	58cm	XL
ST/19	ST/29	52cm	48cm - 56cm	61cm	XXL