

Three Panel Knee Immobiliser

A three panel splint made in high quality foam laminate with bound edges, designed for longer term use. The two detachable side panels can be positioned to fit most leg circumferences. The four steels are removable for the patient to launder the splint.

- + Conservative treatment for knee injuries
- + Maintain knee extension following surgery
- + Maintain knee extension following trauma
- + Maintain knee extension following soft tissue injuries



Code	Splint Length
TPK/12	30cm
TPK/14	36cm
TPK/16	41cm
TPK/18	46cm
TPK/20	51cm
TPK/22	56cm
TPK/24	61cm
TPK/26	66cm
TPK/28	71cm