

# Universal Plantar Faciitis Orthosis

Holds foot in dorsiflexion and allows load application to soft tissues of plantar fascia and achilles tendon. Dorsiflexion easy to adjust using double hook/loop straps additional wedges (3 supplied).

Fits left and right. Replacement liners available.



Code	Replacement Liner Code	Shoe Size	Size
UFO/S	UFO/S-L	3 - 5.5	S
UFO/M	UFO/M-L	6 - 9	M
UFO/L	UFO/L-L	9.5 - 11	L