

Ventilated SK Knee Brace

Using multi-perforated towel-lined latex free neoprene - this knee brace provides excellent medial and lateral support and great patient comfort. The ventilated material reduces moisture build up whilst the towel lining improves comfort

The brace and hinge casings are fully bound to enhance the overall look of the brace whilst the hook/loop straps have grip-taps to assist patients with limited dexterity.

- + Mild instabilities of the knee
- + Meniscus injuries
- + Trauma or degenerative deformities of the knee joint
- + Protection in the later stages of rehab



Code	Knee Circumference	Size
VSKP/02	25cm - 28cm	XS
VSKP/03	28cm - 34cm	S
VSKP/05	34cm - 40cm	M
VSKP/07	40cm - 46cm	L
VSKP/08	46cm - 56cm	XL
VSKP/09	56cm - 61cm	XXL