

Wristlet

A 15cm wrist support made in Neoprene which retains heat to soothe painful joints. The Wristlet conforms well giving light compression, and can be cut to fit swollen joints without fraying.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Instability
- + Post fracture management





Made from Neoprene material for ultimate comfort, which retains heat to soothe painful joints.



Designed specifically for people with potentially limited mobility, this wrist and thumb brace is easy to put on and take off, even one handed.



The Wristlet can be cut to fit swollen joints without fraying.

Code (Left)	Code (Right)	Size
WL/14	WL/24	S - M
WL/16	WL/26	M - L
WL/LXL/L	WL/LXL/R	L - XL