

Wiltshire Maternity Belt

A pelvic support belt with an elastic back section and two padded side straps, plus two variable tension straps.

The slimline design of this belt makes it comfortable to wear when sitting. Provides strong circumferential support around the pelvis holding the two sacroiliac joints and the symphysis pubis joint firmly.

Back depth 10cm, front depth 8cm.

Also available in black.



Code (White)	Code (Black)	Circumference (at widest part of abdomen)
WMB/1	WMB/1/BK	96cm - 116cm
WMB/2	WMB/2/BK	116cm - 130cm
WMB/3	WMB/3/BK	130cm - 146cm
WMB/4	WMB/4/BK	146cm - 160cm