

# Elastic Wrist Thumb Brace

The Elastic Wrist Thumb Brace holds the thumb firmly in abduction, supporting the thumb MP joint and restricting movement at the CMC joint, which is crucial for alleviating pain and assisting the healing process. This thumb support also incorporates a padded aluminium thumb stay and is made from a gently elasticated material that provides compression to the injured area, reducing swelling.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Osteoarthritis
- + Carpal tunnel syndrome
- + CMC joint instability
- + De Quervain's syndrome
- + Gamekeeper's thumb
- + Guyon's canal syndrome
- + Inflamed tendons
- + Luxation
- + Medial subluxation
- + Repetitive strain injury + more






Made from Neoprene material for ultimate comfort, while remaining highly supportive.



Incorporates a padded aluminium thumb stay that provides excellent joint stabilisation. Holds thumb in abduction, supporting MP joint and restricting CMC joint movement.



Designed specifically for people with potentially limited mobility, this wrist and thumb brace is easy to put on and take off, even one handed.

Code (Left)	Code (Right)	Size
WTE/12	WTE/22	XS
WTE/13	WTE/23	S
WTE/15	WTE/25	M
WTE/17	WTE/27	L
WTE/18	WTE/28	XL