

Neoprene Wrist Thumb Brace


Features a malleable aluminium palmar bar that provides fantastic support and stabilisation for the wrist and thumb. This brace also features a specially formulated Neoprene fabric that is ideal for orthotic braces because it provides the perfect combination of insulation and compression offering relief for arthritic conditions and reducing swelling by improving circulation.

This brace both supports the thumb MCP joint and restricts movement at the CMC joint, allowing you to go about your daily life in comfort and reduced pain.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + CMC joint arthritis
- + Carpal tunnel syndrome
- + De Quervain's syndrome
- + Degenerative irritation
- + Inflamed tendons
- + Soft tissue injury
- + Wrist trauma





This brace both supports the thumb MCP joint and restricts movement at the CMC joint, allowing you to go about your daily life in comfort and reduced pain.



Neoprene is an excellent material for both compression and heat retention, making it perfect for use on supports and splints. The neoprene of this support will retain body heat in the wrist, keeping the area warm and increasing the blood circulation to promote proper healing.



| Code (Left) | Code (Right) | Size |
|-------------|--------------|------|
| WTO/12 | WTO/22 | XS |
| WTO/13 | WTO/23 | S |
| WTO/15 | WTO/25 | M |
| WTO/17 | WTO/27 | L |
| WTO/18 | WTO/28 | XL |
| WTO/19 | WTO/29 | XXL |