

# Air X


## Two Piece Wrist

This two piece brace is shaped to achieve excellent fit, compression and stability around the hand, wrist and forearm. The 3D structure of the AirX textile gives a high degree of cushioning and increases comfort and compression. It also allows air transfer making the brace comfortable to wear for extended periods.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Instability
- + Post fracture management





Air X material is soft, lightweight and perfectly follows the body's contours. The moisture-transpiring fabric prevents sweat build-up and keeps skin feeling comfortable.



Reinforced bar casing prevents any chafing and adds additional stabilisation, with a rigid aluminium palmar bar to provide additional support.

Code (Left)	Code (Right)	Size
AFS/12	AFS/22	XS
AFS/13	AFS/23	S
AFS/15	AFS/25	M
AFS/17	AFS/27	L
AFS/18	AFS/28	XL