

Ascender Telescopic ROM Elbow Orthosis

The Ascender Telescopic ROM Elbow Orthosis is easily adjustable and secure. With additional features offering superior comfort and support to patients recovering from surgery, injuries or instabilities.

The range of motion settings can be easily adjusted via the spring loaded quick-clips by a simple 'pull and slide' motion. The design prevents unintentional changes to the flexion and extension settings. The drop lock is ergonomically designed to be easily locked and unlocked due to its curved concave profile. Can be used for an array of injuries and conditions including;

- + Following surgery of the elbow joint.
- + For immobilisation after elbow joint injuries and subsequent mobilisation.
- + Elbow ligament injuries.



Code (Left)	Code (Right)	Code (Left + Supination Bar)	Code (Right + Supination Bar)	Size
ASCE/L/REG	ASCE/L/REG	ASCE/L/REG/SUP	ASCE/R/REG/SUP	Universal