

Lycrafleece Thumb Spica

Supports the wrist and thumb, leaving the palm completely free. An 11cm removable flexible stay restricts movement at the CMC and MCP joints, leaving the IP free. The angled hook/loop attachment at the wrist and thumb ensures a perfect fit. Use for patients with Rheumatoid and Osteoarthritis.

Lycrafleece is a Lycra Spandex 3 ply laminate, designed for the ultimate in patient comfort. The unique construction makes this material perfect for orthopaedic bracing. It is our most comfortable and durable material which makes this range of products the elite in medical bracing.

Lycrafleece has unique properties which keep the patient warm when its cold and cool when its warm for optimum comfort without sacrificing performance.



The angled hook/loop attachment at the wrist and thumb ensures a perfect fit.

| Code (Left) | Code (Right) | Size |
|-------------|--------------|------|
| BLTS/12 | BLTS/22 | XS |
| BLTS/13 | BLTS/23 | S |
| BLTS/15 | BLTS/25 | M |
| BLTS/17 | BLTS/27 | L |
| BLTS/18 | BLTS/28 | XL |