

# Lycrafleece Wrist Brace

An 18cm three-piece wrist brace featuring longitudinal seams to give greater conformity at the wrist. The malleable palmar bar has a spoon shaped end which conforms to the palm of the hand for a comfortable fit.

Lycrafleece is a Lycra Spandex 3 ply laminate, designed for the ultimate in patient comfort. The unique construction makes this material perfect for orthopaedic bracing. It is our most comfortable and durable material which makes this range of products the elite in medical bracing.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Instability
- + Post fracture management





Lycrafleece has unique properties which keep the patient warm when its cold and cool when its warm for optimum comfort without sacrificing performance.



Lycra Spandex gives the material unique stretch and recovery properties whilst a polyester lining material wicks moisture away from the skin to allow greater patient compliance.



The malleable palmar bar has a spoon shaped end which conforms to the palm of the hand for a comfortable fit.

Code (Left)	Code (Right)	Size
BLWB/12	BLWB/22	XS
BLWB/13	BLWB/23	S
BLWB/15	BLWB/25	M
BLWB/17	BLWB/27	ML
BLWB/18	BLWB/28	XL