

Beta Wrist Brace

An excellent value and best seller wrist brace made of a dual compression elastic which conforms to the wrist.

The vinyl pocket encloses a removable aluminium palmar bar which can be shaped to the wrist for a comfortable fit. Lycra binding finishes the edges for comfort and prevents chafing at the thumb web. Available in beige or black.

Provides wrist support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions.
- + Instability
- + Post fracture management



Beige - Left	Beige - Right	Black - Left	Black - Right	Size	Wrist Circ.
BWB/12	BWB/22	BWB/12/BK	BWB/22/BK	XS	> 15cm
BWB/13	BWB/23	BWB/13/BK	BWB/23/BK	S	15cm - 17cm
BWB/15	BWB/25	BWB/15/BK	BWB/25/BK	M	17cm - 19cm
BWB/17	BWB/27	BWB/17/BK	BWB/27/BK	L	19cm - 22cm
BWB/18	BWB/28	BWB/18/BK	BWB/28/BK	XL	23cm - 25cm