

Deltaform Wrist Brace

The Deltaform Wrist Brace offers support, compression and stabilisation. It's made from a special elasticated cotton and features a rigid aluminium bar for extra stabilisation and support. The gentle compression provided by the elastic helps to reduce swelling and assists with the natural healing process.

The simple hook loop fastening, which fastens securely to provide a personalised fit and easy application.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthralgia
- + Arthritis
- + Carpal tunnel syndrome
- + Degenerative diseases
- + Degenerative irritation
- + General weakness/instability
- + Guyon's canal syndrome
- + Inflamed/swollen tendons
- + Medial subluxation
- + Overuse injuries
- + Post carpal/metacarpal surgery
- + Post surgical repair
- + Soft tissue injury
- + Sprained wrist
- + Tendinosis
- + Wrist trauma





Comfortable elastic/cotton material and longitudinal seam for excellent conformity to the wrist and forearm.



Comes in a choice of beige or black colours with durable protection is ideal for longevity.



Features a removable metal stay, which provides extra support and is malleable and can be adjusted to provide the perfect fit.

Code (Left Black)	Code (Right Black)	Code (Left Beige)	Code (Right Beige)	Size
DFS/12/BK	DFS/22/BK	DFS/12	DFS/22	XS
DFS/13/BK	DFS/23/BK	DFS/13	DFS/23	S
DFS/15/BK	DFS/25/BK	DFS/15	DFS/25	M
DFS/17/BK	DFS/27/BK	DFS/17	DFS/27	L
DFS/18/BK	DFS/28/BK	DFS/18	DFS/28	XL