

Freedom Wrist Brace

An 18cm long brace with cut-away design to allow flexion at the palmar crease and opposition of the thumb to the fingers. The soft Lycra binding prevents chafing of the thumb web. Longitudinal seam gives good conformity and support to the wrist and forearm.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Soft tissue injuries and conditions
- + Instability
- + Post fracture management



Code (Left Beige)	Code (Right Beige)	Code (Left Black)	Code (Right Black)	Size
FWB/12	FWB/22	FWB/12/BK	FWB/22/BK	XS
FWB/13	FWB/23	FWB/13/BK	FWB/23/BK	S
FWB/15	FWB/25	FWB/15/BK	FWB/25/BK	M
FWB/17	FWB/27	FWB/17/BK	FWB/27/BK	L
FWB/18	FWB/28	FWB/18/BK	FWB/28/BK	XL



The soft Lycra binding prevents chafing of the thumb web and a longitudinal seam to give good conformity.



Reinforced bar casing prevents any chafing and adds additional stabilisation, with a rigid aluminium palmar bar to provide additional support.

Code (Left Beige)	Code (Right Beige)	Code (Left Black)	Code (Right Black)	Size
FWB/12	FWB/22	FWB/12/BK	FWB/22/BK	XS
FWB/13	FWB/23	FWB/13/BK	FWB/23/BK	S
FWB/15	FWB/25	FWB/15/BK	FWB/25/BK	M
FWB/17	FWB/27	FWB/17/BK	FWB/27/BK	L
FWB/18	FWB/28	FWB/18/BK	FWB/28/BK	XL