

FixStep Walker

NEW PRODUCT

Our new and improved fixed walker available in standard and short sizes with high volume pump air liners or standard liners.

- + Excellent rocker sole profile for improved gait
- + Two-part injection moulded sole unit
- + Improved shock absorbing - less impact during walking
- + Flared & deeper heel section to retain foot and reduce any pressure on Achilles area
- + High quality hook/loop straps - easy for patients to open
- + Strap lock system avoids straps rotating
- + Moulded high volume pump on air liner for faster inflation and easier patient use
- + Comfortable detachable liners with replacements available
- + Reinforced moulded plastic uprights - metal free - making them x-ray safe



Size	FixStep with Standard Liner	Replacement Standard Liner for FixStep	Short FixStep with Standard Liner	Replacement Standard Liner for Short FixStep	FixStep with Air Liner	Short FixStep with Air Liner	Shoe Size (UK)
XS	FXP/02	FXPL/02	FXPS/02	FXPSL/02	FXPA/02	FXPSA/02	1.5 - 3
S	FXP/03	FXPL/03	FXPS/03	FXPSL/03	FXPA/03	FXPSA/03	3 - 5.5
M	FXP/05	FXPL/05	FXPS/05	FXPSL/05	FXPA/05	FXPSA/05	6 - 9
L	FXP/07	FXPL/07	FXPS/07	FXPSL/07	FXPA/07	FXPSA/07	9.5 - 11
XL	FXP/08	FXPL/08	FXPS/08	FXPSL/08	FXPA/08	FXPSA/08	11.5 - 13.5

LITFXPSELB 12/23