

Juraprene Long Wrist Wrap

A latex free neoprene with terry plush lining makes this a very comfortable support. Full finger and thumb flexion whilst wider reinforced hook and loop closure gives improved circumferential compression and support to the wrist. Available in Beige, Black and Blue.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Sprains and strains
- + Instability
- + Post fracture management





Neoprene fabric is gently elastic, breathable and durable. The terry plush lining provided makes the inside of the wrist wrap very comfortable.



Full finger and thumb flexion means minimal interference with your day-to-day routine, plus the reinforced hook and loop closures make this wrist wrap simple and secure to fasten.

Code (Beige)	Code (Blue)	Code (Black)	Size
JLWW/03/BE	JLWW/03/BL	JLWW/03/BK	S
JLWW/05/BE	JLWW/05/BL	JLWW/05/BK	M
JLWW/07/BE	JLWW/07/BL	JLWW/07/BK	L