

# Jura Fixed Walker

One piece injection moulded walker with 'paddle' uprights to reduce pressure on lower limb and provide improved stability and immobilisation. Foot section can accommodate full dressings and padding.

New lightweight design sole unit to greatly improve shock absorption, gait pattern to reduce the overall weight of the walker to an absolute minimum.

Reduced strap fastening to simplify removal/application and improve patient acceptance. Comfortable detachable liners with replacements available.

- + Foot and ankle conditions that benefit from motion restriction
- + Ankle sprains
- + Trauma and rehabilitation
- + Stable foot fractures



Code (Tall)	Code (Short)	Code (Air Liner - Tall)	Code (Air Liner - Short)	Shoe Size (UK)	Size
JWF/02	JWFS/02	JWF/02/A	JWFS/02/A	1.5 - 3	XS
JWF/03	JWFS/03	JWF/03/A	JWFS/03/A	3 - 5.5	S
JWF/05	JWFS/05	JWF/05/A	JWFS/05/A	6 - 9	M
JWF/07	JWFS/07	JWF/07/A	JWFS/07/A	9.5 - 11	L
JWF/08	JWFS/08	JWF/08/A	JWFS/08/A	11.5 - 13.5	XL