

# Jura Wrist / Thumb Brace Long

A wrist thumb orthosis with malleable aluminium thumb stay to hold the thumb firmly in abduction, the Jura Wrist / Thumb Brace supports the thumb MP joint and restricts movement at the CMC joint.

Provides wrist support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis.
- + Sprains and strains.
- + Soft tissue injuries and conditions.
- + Instability.
- + Post fracture management.



Code (Left)	Code (Right)	Size	Wrist Circumference
JWTL/12	JWTL/22	XS	< 15cm
JWTL/13	JWTL/23	S	15 - 17cm
JWTL/15	JWTL/25	M	17 - 19cm
JWTL/17	JWTL/27	L	19 - 22cm
JWTL/18	JWTL/28	XL	23 - 25cm