

Nuform Wrist Brace

An 18cm wrist brace made from ventilated elastic, which provides gentle compression and excellent breathability. This wrist brace features a vinyl pocket with an aluminium palmar bar for increased support and stabilisation, which can be shaped to provide a personalised fit. Available in beige or black.

This orthopaedic support offers a fantastic stabilisation at the wrist. By stabilising the joint, the wrist brace ensures that there is no excessive movement relative to the surrounding structures.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Arthritis
- + Carpal tunnel syndrome
- + Inflamed/swollen tendons
- + Soft tissue injury
- + Degenerative irritation
- + Tendonitis
- + Worn cartilage
- + Wrist trauma





The support offered by this brace combined with the unrestricted finger and thumb movement allows you to complete daily tasks without interrupting the healing process.



Comfortable and ventilated elastic provides gentle compression and excellent breathability.

Code (Left)	Code (Right)	Size
NUF/12	NUF/22	XS
NUF/13	NUF/23	S
NUF/15	NUF/25	M
NUF/17	NUF/27	L
NUF/18	NUF/28	XL