


Off-Loading Insole

The off-loading insole has a self-adhesive stabilising base board, with removable sections of the insole for custom off-loading in treatment of various foot pathologies.

The insoles greatly reduces pressure on plantar surface when patient is walking. Available for all Jura post-op shoes.



Code (Mens)	Shoe Size (Mens)	Code (Womens)	Shoe Size (Womens)	Size
P-OP/MSI	5 - 7	P-OP/WSI	3 - 5	S
P-OP/MMI	7.5 - 9	P-OP/WMI	5.5 - 7	M
P-OP/MLI	9.5 - 11	P-OP/WLI	7.5 - 9	L
P-OP/XLI	11.5 - 13	-	-	XL