

Jura Post-Op Shoe

The square toe configuration of the Jura Post-Op Shoe offers protection and improved fit/ accommodation. The rocker-sole profile aids gait and reduces forces on foot. Outsole is non-slip and material composition provides reduced movement under the metatarsal heads. Simple closure system allows for additional dressings. Insole can be removed for alteration, improved accommodation or fitting 'Off-Loading Insole'. Toe cover can be removed if required.



- + Accommodate bandages or swelling
- + Foot surgery
- + Broken toe management
- + Metatarsal injuries
- + Soft tissue procedures
- + Management of diabetic foot
- + Post-operative recovery



Code (Mens)	Shoe Size (Mens)	Code (Womens)	Shoe Size (Womens)	Size
P-OP/MS	5 - 7	P-OP/WS	3 - 5	S
P-OP/MM	7.5 - 9	P-OP/WM	5.5 - 7	M
P-OP/ML	9.5 - 11	P-OP/WL	7.5 - 9	L
P-OP/XL	11.5 - 13	-	-	XL

Removable Toe Cover Code	Size
P-OP/MS	Universal