

Procool D-Ring Wrist Splint

Low profile support with circumferential slip-on design and D-Ring fastening to allow easier one handed application. Aluminium palmar bar can be removed for further adjustment. Dorsal stockinette lining facilitates application and cushions strap closures.

Design of the Standard and Firm braces allow full finger and thumb function whilst limiting wrist flexion.

The Firm splint has the addition of adjustable dorsal stays to reduce unwanted flexion and extension and improve overall immobilisation of acute soft tissue injuries. This brace can be trimmed if required.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Sprains
- + Arthritis
- + Carpal Tunnel Syndrome
- + Strains





Slip-on design and D-ring fastening to allow easier one handed application.



The design of the Standard and Firm braces allow full finger and thumb function whilst limiting wrist flexion.



Reinforced bar casing prevents any chafing and adds additional stabilisation, with a rigid aluminium palmar bar to provide additional support.

Code (Left)	Code (Right)	Code (Left Firm)	Code (Right Firm)	Size	Wrist Circumference
PDR/12	PDR/22	PDRF/12	PDRF/22	XS	> 15cm
PDR/13	PDR/23	PDRF/13	PDRF/23	S	15 - 17cm
PDR/15	PDR/25	PDRF/15	PDRF/25	M	17 - 19cm
PDR/17	PDR/27	PDRF/17	PDRF/27	L	19 - 22cm
PDR/18	PDR/28	PDRF/18	PDRF/28	XL	23 - 25cm