

Procool Deluxe Wrist Thumb Brace

A uniquely shaped wrist thumb brace offering excellent support and compression to the wrist, CMC and MCP joints. Malleable palmar bar positions wrist in functional extension and flexible stay provides abduction and limited motion to support and protect CMC & MCP joints. Adjustable thumb closure allows easy application of the brace.

Procool is a range of lightweight and breathable splints producing excellent compression and support. Procool offers you just the fit you want - excellent stability without excessive bulk.

Provides wrist and thumb support for an array of injuries, conditions, pre-and post-operations such as:

- + Carpal Tunnel Syndrome
- + De Quervain's
- + Post-scaploid Fracture
- + Post-surgical Recovery
- + Repetitive Strain Injury
- + Skier's Thumb
- + Tendinitis



Incorporates a padded aluminium thumb stay for great protection and support, and the aluminium bar through the palm keeps the wrist stable.



Made with tiny perforations in the Neoprene for comfortable wear and a terry cloth lining that feels good against the skin.



Adjustable thumb closure allows for easy application, protecting and supporting without restricting necessary movement. The brace also comes complete with three thick hook / loop straps which allow the user to customise tightness based on the requirements of their condition.



The brace comes complete with three thick hook / loop straps which allow the user to customise tightness based on the requirements of their condition.

Code (Left)	Code (Right)	Size
PDWT/12	PDWT/22	XS
PDWT/13	PDWT/23	S
PDWT/15	PDWT/25	M
PDWT/17	PDWT/27	L
PDWT/18	PDWT/28	XL